

Gabion Tensile Grips

Introduction

- Used on electronic universal testing machine for determination of tensile strength for wire mesh and selvedge strength for wire mesh at room temperature.
- Principle: The specimen is stretched using a tensile testing machine until the mesh surface breaks or the mesh edge flares open or the steel wire breaks, in order to determine the mechanical properties of the mesh surface.



For further information please contact:

T_{EST}
M_{ACHINES}
A_{USTRALIA}

0418 369 505

sales@testmachines.com.au

www.testmachines.com.au